

# DS2 ●

## One of a kind

The **DS-2** is a **unique Digital monoscope**, based on the conventional model MS-2. With its long working distance it is perfect for viewing macro specimens such as plants, insects and minerals.

The **LED illumination**, mounted on a flexible gooseneck arm, guarantees sufficient light from different angles. The **LED illumination** works with rechargeable batteries and can therefore also be used during outdoor biology classes or studies in the backyard.

By connecting the **USB 2.0** cable to your computer and using one of the software provided with the microscope, **Motic Play** or **Motic Educator**, you can **easily capture digital images or create video's**. Besides the Motic Play and Educator, we also offer the possibility to download the sophisticated **Motic Images Plus 2.0** software from our website.

## GOINGgreen



- Monocular head with built-in digital camera CMOS 1/4" chip
- 800x600 Pixels USB2.0
- Widefield eyepiece WF10X/18mm
- 2X built-in monoscope objective + additional 2x attachable
- Pole-Type stand with slow falling focusing mount
- Reversible stage plate black/white
- Rechargeable LED illumination mounted on flexible gooseneck
- Safety auto cut-off charger and batteries included
- Main supply 100V-240V (CE)
- Motic Play for PC & Mac
- Motic Educator for PC & Mac
- Motic Images Plus 2.0 as download
- Calibration slide and dust cover

